# GSD NovaType SARS-CoV-2 ID

## RT-PCR

**RUO**

**FOR RESEARCH USE ONLY**  
**Not for use in diagnostic procedures**

English .....	2
Deutsch .....	6
Français .....	11
Italiano .....	16
Español .....	21
Português .....	26
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía / Bibliografia .....	31
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas .....	31
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos .....	32

---

Product Number:	PCOV6050T	(1 x 96 Determinations)
	PCOV6054T	(4 x 96 Determinations)

---

## ENGLISH

### 1. INTRODUCTION

---

End of 2019, a novel respiratory disease emerged in the city of Wuhan, Hubei Province of the People's Republic of China, and soon spread rapidly within the country and worldwide. The causative agent was identified as severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). SARS-CoV-2 (2019-nCoV), like the closely related SARS coronavirus (SARS-CoV), belongs to the genus Betacoronavirus within the family of coronaviruses. The zoonotic reservoir of the virus appears to be bats.

Coronaviruses are enveloped, positive single-stranded large RNA viruses that infect humans, but also a wide range of animals. The common human coronaviruses NL63, 229E, OC43 and HKU1 are widespread especially throughout the winter months. They are responsible for up to one third of all acute respiratory diseases, typically with mild symptoms (common cold). More than 80 % of the adult population have antibodies against human coronaviruses. The immunity from previous infections lasts only for a short period of time. Therefore, reinfections with the same pathogen are possible just after one year.

SARS-CoV-2 is predominantly transmitted by droplet infection via coughing or sneezing and through close contact with infected patients. In theory, smear infection and infection through the conjunctiva of the eyes are also possible. The incubation period is in the median 5–6 days (and up to 14 days maximum).

The clinical manifestations of SARS-CoV-2-related COVID-19 disease include fever, cough, respiratory problems and fatigue. In most patients the infection manifests with symptoms of a mild febrile illness with irregular lung infiltrates.

The initial clinical sign of COVID-19 which allowed case detection was pneumonia. But it turned out that the course of the disease is non-specific and varies widely, from asymptomatic courses to severe pneumonia with lung failure and death. However, based on current knowledge, around 80 % of the illnesses are mild to moderate.

Although severe courses of the disease also occur in younger patients and people without previous illness, the following groups of people have an increased risk of serious forms of the disease: elderly people (with a steadily increasing risk from around 50-60 years of age), smokers and people with certain diseases of the cardiovascular system or the lungs, patients with chronic liver diseases, diabetes mellitus, cancer, or patients with a weakened immune system (e.g. due to immune deficiencies or by taking drugs that suppress the immune system).

Treatment focuses on supportive measures depending on the severity of the clinical picture (e.g. oxygen administration, fluid balance management, etc.) as well as the treatment of relevant underlying diseases. Various specific therapeutic approaches (directly antiviral effective, immunomodulatory effective) have been and are being investigated in studies during the course of the SARS-CoV-2 pandemic.

In late 2020 two COVID-19 variants were identified that show mutations in different positions of the S gene, coding for the spike protein. Some of these mutations result in higher infectivity. In particular mutation N501Y is identified as one root cause. This mutation is found in virus isolates in United Kingdom (UK), lineage B.1.1.7, and in South Africa, lineage B.1.351. These two variants can be distinguished by mutation A570D which is additionally present only in the variant from UK.

Species	Disease	Symptoms e.g.	Transmission route
SARS-CoV-2 (severe acute respiratory syndrome coronavirus 2)	COVID-19	the course of the disease is unspecific, diverse and varies greatly, from asymptomatic courses to severe pneumonia with lung failure and death	primary mode of transmission: droplet infection; smear infections and infections via the conjunctiva of the eyes are theoretically possible

Infection or presence of pathogen may be identified by

- Nucleic acid Amplification Techniques (NAT): e.g. RT-PCR
- Serology: detection of antibodies by e.g. ELISA

### 2. INTENDED USE (FOR RESEARCH USE ONLY!)

---

The GSD NovaType SARS-CoV-2 ID is a research use only (RUO) real-time reverse transcription PCR assay designed for the simultaneous qualitative detection and discrimination of the SARS-CoV-2 (Severe acute respiratory syndrome coronavirus 2) S gene variants wildtype, and mutations N501Y and A570D. Sample material should be genomic RNA extracted from human respiratory (nasal wash/smear, nasopharyngeal wash/smear, oropharyngeal swab, bronchoalveolar lavage, and gargle solution) specimen types that have previously tested positive for SARS-CoV-2 RNA by diagnostic RT-PCR methods e.g. GSD NovaPrime® SARS-CoV-2 (COVID-19). The GSD NovaType SARS-CoV-2 ID is expressly not intended for diagnosis of SARS-CoV-2 infection and not for use in diagnostic procedures. It is only designed to aid investigations related to the prevalence and frequency of virus variants.

### 3. PRINCIPLE OF THE ASSAY

---

The qualitative determination of specific RNA is based on Real-Time reverse-transcription Polymerase Chain Reaction (RT-PCR) technology. The kit contains specific primers and probes labelled with fluorescent reporter and quencher dyes for amplification and simultaneous detection of specific RNA sequences which represent **specific SARS-CoV-2 S gene variants**.

The gene of interest specific probes are labelled with the fluorophores HEX™ and FAM™ thereby allowing parallel discrimination between S gene wildtype and mutations N501Y and A570D.

The GSD NovaType SARS-CoV-2 ID RT-PCR has been validated for following real-time PCR platforms:

- ABI Prism® 7500 SDS (Applied Biosystems™)
- ABI Prism® 7500 Fast SDS (Applied Biosystems™)
- AriaMx™ (Agilent Technologies)
- AriaDx™ (Agilent Technologies)

## 4. MATERIALS

### 4.1. Reagents Provided

Cap	Symbol	Component	PCOV6050T (PCOV6054T)	
			Number of Vials	Volume per Vial [µL]
green	E-MIX	RT-PCR Enzyme Mix (reverse transcriptase, hot-start DNA polymerase, RNase inhibitor, nucleotides, magnesium, stabilizers and buffer)	1 (4)	550
blue/white	PP VAR	Primer-Probe-Mix (for discrimination of variants)	1 (4)	350

### 4.2. Materials and Equipment needed, but not provided

- Deionized nuclease-free and nucleic-acid free water
- Real-Time PCR instrument (for already validated instruments refer to 3.PRINCIPLE OF THE ASSAY)  
Alternative Real-Time PCR instruments might also be appropriate. Their suitability for use with the GSD NovaType SARS-CoV-2 ID has to be validated by the user.
- Appropriate Real-Time PCR consumables (e.g. disposable tubes, reaction plates, corresponding optical closing materials)
- Benchtop microcentrifuge
- Centrifuge with a rotor for microtiter plates
- Vortex mixer
- Adjustable pipettes in relation to reaction setup
- Disposable DNase/RNase free pipette tips with filters
- Disposable powder-free gloves

## 5. STABILITY AND STORAGE

The GSD NovaType SARS-CoV-2 ID kit is shipped on dry ice and all components should arrive frozen.

- All components have to be stored between -30 °C and -15 °C immediately after arrival.
- Repeated freeze thaw cycles (more than two) of reagents should be avoided, since this might affect the performance of the kit. Reagents should be frozen in aliquots if they are used intermittently.
- Keep unfrozen storage (e.g. storage on ice) as short as possible.
- Keep the E-MIX and the PP|VAR in the freezer, until you are ready to use it.
- Protect the E-MIX and the PP|VAR from light.

## 6. SAMPLE PREPARATION

- Extracted RNA or total nucleic acid extracted from human respiratory specimen types (nasal lavage/smear, nasopharyngeal lavage/smear, oropharyngeal swab, bronchoalveolar lavage, and gargle solution) previously tested positive for SARS-CoV-2 RNA by RT-PCR methods e.g. GSD NovaPrime® SARS-CoV-2 (COVID-19) is the starting material for the GSD NovaType SARS-CoV-2 ID.
- The quality of the extracted RNA has a crucial effect on the performance of the entire RT-PCR test system. Make sure that the nucleic acid extraction method is compatible with Real-Time PCR technology.

## 7. ASSAY PROCEDURE

### 7.1. Reaction Setup

- Please read the instructions for use carefully **before** performing the assay. Reliability of results depends on following strictly the instructions for use.
- Before use make sure that all samples and reagents are thawed completely, mixed by up and down pipetting or vortexing and centrifuged briefly.
- Important: E-MIX is viscous. Briefly spin down the tube to ensure material has not lodged in the cap or side of tube. Be sure to pipette and dispense carefully, and use pipette tips suitable for pipetting viscous liquids.
- The use of nuclease-free and nucleic-acid free water as no template control (NTC) is highly recommended.
- Define the positions (wells) for samples and controls (NTC) on the plate.

Reaction Setup	
E-MIX	5 µL
PP VAR	3 µL
Sample or NTC	12 µL
<b>Total volume</b>	<b>20 µL</b>

- Close the optical reaction plate with corresponding optical closing material.
- Centrifuge the optical reaction plate in a centrifuge with a rotor for microtiter plates for 60 seconds at approximately 1000 x g (~ 3000 rpm).

## 7.2. Programming the Real-Time PCR Instrument

Regarding setup and programming of the Real-Time PCR instrument, please use the manual of the respective instrument.

For detailed programming instructions regarding the use of the GSD NovaType SARS-CoV-2 ID RT-PCR on specific Real-Time PCR instruments please contact NovaTec Immundiagnostica GmbH.

RT-PCR Run Settings	
Reaction Volume	20 µL
Ramp Rate	Standard or Fast
Passive Reference	ROX™

### Fluorescent Detectors / Fluorophores

Detection of the amplified viral nucleic acid fragment is performed in the detection channels FAM™ (S gene mutation N501Y and A570D) and HEX™ (S gene wildtype N501Y and A570D) with MGBEQ Quencher. Settings for the ABI Prism® 7500 SDS/Fast SDS (Applied Biosystems™) are given in brackets.

Detection	S Gene Variant	Fluorophore	Quencher
SARS-CoV-2 S gene	N501Y mutation	FAM™	MGBEQ (none)
	N501 wildtype	HEX™ (VIC®)	
	A570D mutation	FAM™	
	A570 wildtype	HEX™ (VIC®)	

### Temperature Profile and Data Collection

No. of Cycles	Temperature	Time (min)	Data Collection
1	50 °C	10:00	-
1	95 °C	03:00	-
40	95 °C	00:10	-
	60 °C	00:30	Fluorescence measurement at the end of every cycle

Before starting the test run, please check the settings for cycles, temperature and time.

## 8. RESULTS

Data analysis should be performed with the software of the used real-time PCR device according to manufacturer's instructions.

Analysis settings:

Setting	Recommendation
Threshold	Enter a value for the threshold so that the threshold is: <ul style="list-style-type: none"> <li>• Above the background</li> <li>• Below the plateau and linear regions of the amplification curve</li> <li>• Within the exponential phase of the amplification curve</li> </ul>
Baseline	Select start and end cycle values so that the baseline ends before significant fluorescence is detected.

### 8.1. Run Validation Criteria

Test run is valid only if RT-PCR run complete.

The NTC Ct values must meet the acceptance criterion in the table below for the assay to be valid.

Only after that sample data may be analyzed and interpreted.

Validation Criteria	Result/Acceptable Ct	Valid/Invalid	Measure
NTC	not detected (ND)	valid	-
	Ct ≤ 40	invalid	Repeat test run

## 8.2. Interpretation of Results

If test run is valid interpretation of sample results is as follows:

- Positive (POS): Ct ≤ 38
- Negative (neg): Not detected (ND) or Ct > 38

SARS-CoV-2 S gene	Detection Channel	
	HEX™ (VIC®)	FAM™
wildtype	POS	neg
lineage B.1.1.7 (N501Y/A570D)	neg	POS
lineage B.1.351 (N501Y)	POS	POS

Discrimination of virus variants with respect to S gene codons 501 and 570 is only possible for samples

- that have already been pretested positive for the presence of SARS-CoV-2 RNA and
- whose Ct value in one or both of the S gene channels is below Ct ≤ 38

## 9. LIMITATIONS OF THE PROCEDURE

Very low amounts of starting material, e.g. weak positive samples (high Ct values in the previous diagnostic SARS-CoV-2 PCR), can result in GSD NovaType SARS-CoV-2 ID test runs with negative signals. However, such results do not exclude the presence of SARS-CoV-2 RNA.

## 10. QUALITY CONTROL

In accordance with NovaTec Immundiagnostica GmbH ISO-certified Quality Management System, each lot of GSD NovaType SARS-CoV-2 ID has been tested against predetermined specifications to ensure consistent product quality.

## 11. TRADEMARKS AND DISCLAIMERS

Registered names, trademarks, etc. used in this document are to be considered protected by law even if not specifically marked as such.

## 12. PRECAUTIONS AND WARNINGS

- **FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.**
- All human samples should be regarded and handled as potentially infectious.
- Do not interchange reagents of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Wear disposable powder-free gloves, a laboratory coat and eye protection when handling specimens.
- Always use DNase/RNase-free disposable reaction tubes and pipette tips with aerosol barriers.
- Avoid microbial and nuclease (DNase/RNase) contamination of the specimen and the components of the kit.
- In order to avoid contamination of working space with nucleic acids, reaction tubes/plates should not be opened after amplification.
- RT-PCR is highly sensitive to nucleic acid contamination. Therefore, positive / potentially positive material must be stored separately from all other components of the kit.
- Dedicate supplies and equipment to the separate working areas and do not move them from one area to another.
- This assay must not be used on the specimen directly.
- Prior to using this assay the nucleic acid has to be extracted with suitable extraction methods from the original specimen.
- Since ethanol is a strong Real-Time PCR inhibitor, it is necessary to completely eliminate it prior to the elution of the nucleic acid during extraction. If using spin columns with washing buffers containing ethanol, it is highly recommended to perform an additional centrifugation step of 10 min at approximately 17000 x g (~ 13000 rpm) before eluting the RNA. For this additional centrifugation step, use a new collection tube.
- The result of this RT-PCR kit may be influenced by potential mutations in the genome of the pathogen if they are located in the primer / probe binding region. Underestimation and/or failure to detect the pathogen may occur.
- RT-PCR inhibitors may also elicit underestimation, false negative results or invalid runs. Therefore, only use nucleic acids extraction kits, which remove RT-PCR inhibitors and which are dedicated for downstream RT-PCR processes.
- The Real-Time PCR is only designed for qualified personnel who are familiar with good laboratory practice and trained in RT-PCR.
- For further internal quality control each laboratory should additionally use known samples.

### 12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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